



User's Manual

Porcine Clara cell Specific protein CC16 ELISA Kit



DEIA-BJ2756



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

Porcine Clara cell Specific protein CC16 ELISA Kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of the Clara cell Specific protein CC16. This ELISA kit is for research use only, not for therapeutic or diagnostic applications.

Storage

All components of this kit are stable at 2-8°C until the kit's expiration date.

Detection Range

250-5000 pg/mL

Sensitivity

1.0 pg/mL