



User's Manual

Monkey Defensin Alpha 5, Paneth Cell Specific ELISA Kit



DEIA-BJ2895



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**



Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)



Fax: 1-631-938-8221



Email: info@creative-diagnostics.com



Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

Monkey Defensin Alpha 5, Paneth Cell Specific ELISA Kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of the Defensin Alpha 5, Paneth Cell Specific. This ELISA kit is for research use only, not for therapeutic or diagnostic applications.

Storage

All components of this kit are stable at 2-8°C until the kit's expiration date.

Detection Range

100-2500 pg/mL

Sensitivity

1.0 pg/mL