



User's Manual

Monkey Inhibitor of kB kinase α ELISA Kit

REF DEIA-BJ2900

 96T



RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

Monkey Inhibitor of κ B kinase α ELISA Kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of the Inhibitor of κ B kinase α ELISA. This ELISA kit is for research use only, not for therapeutic or diagnostic applications.

Storage

All components of this kit are stable at 2-8°C until the kit's expiration date.

Detection Range

5-100 ng/mL

Sensitivity

0.1 ng/mL