



User's Manual

Human Phosphorylation P38 Mitogen-activated Protein Kinase ELISA kit

REF

DEIA-BJ344



96T



RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address:** 45-1 Ramsey Road, Shirley, NY 11967, USA

 **Tel:** 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  **Fax:** 1-631-938-8221

 **Email:** info@creative-diagnostics.com  **Web:** www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

Human Phosphorylation P38 Mitogen-activated Protein Kinase ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of the Phosphorylation P38 Mitogen-activated Protein Kinase. This ELISA kit is for research use only, not for therapeutic or diagnostic applications.

Storage

All components of this kit are stable at 2-8°C until the kit's expiration date.

Detection Range

50-1000 ng/mL

Sensitivity

1.0 ng/mL