



User's Manual

Human α Synuclein ELISA kit



DEIA-BJ866



96T





This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

For quantitative detection of SNCA (Alpha-synuclein) in serum, plasma, tissue homogenates and other biological fluids.

Storage

All components of this kit are stable at 2-8°C until the kit's expiration date.

Precision

Intra-Assay: CV<8%

Inter-Assay: CV<10%

Detection Range

15.6-1000 pg/ml

Sensitivity

<9.375 pg/ml

