



User's Manual

Bovine Interleukin 1 ELISA Kit



DEIA-CL012



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

This IL-1 ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Bovine IL-1. This ELISA kit for research use only.

Storage

Store all reagents at 2-8°C.

Sensitivity

1.0 pg/mL