



User's Manual

Bovine Insulin like Growth Factor 1 ELISA Kit



DEIA-CL014



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

This IGF-1 ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Bovine IGF-1. This ELISA kit for research use only.

Storage

Store all reagents at 2-8°C.

Sensitivity

1.0 pg/mL