



User's Manual

Human CVB- IgG(Coxsackievirus B- Immunoglobulin G) ELISA Kit

REF

DEIA-NS2307-102



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

This immunoassay kit allows for the qualitative determination of CVB-IgG in human serum or plasma.

Storage

2-8 °C for 6 months

Precision

Intra-Assay: CV<8%

Inter-Assay: CV<10%