



User's Manual

Human CMV-IgG(cytomegalovirus- Immunoglobulin G) ELISA Kit



DEIA-NS2307-106



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

CMV-IgG ELISA Kit allows for the in vitro quantitative determination of CMV-IgG concentrations in serum, plasma, tissue homogenates and other biological fluids.

Storage

2-8 °C for 6 months

Precision

Intra-Assay: CV<8%

Inter-Assay: CV<10%

Detection Range

Qualitative

Sensitivity

< Qualitative