



User's Manual

Tyrosine Colorimetric Assay Kit



DEIA-PYR064



100T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The Tyrosine Colorimetric Assay Kit is intended for the detection of tyrosine in rat serum, plasma, urine and other body fluids.

Storage

Store kit at -20°C, protected from light. Briefly centrifuge small vials prior to opening.