



## User's Manual

# CCL21 ELISA Kit

REF

DEIA-XY102



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

Enzyme-linked immunosorbent assay (ELISA) for the determination of CCL21 in cell culture supernatants, serum, or plasma.

### Storage

2-8°C

### Detection Range

63 pg/mL-4,000 pg/mL

### Sensitivity

The minimum detectable dose of this CCL21 ELISA typically is less than 12 pg/mL.

### Specificity

This CCL21 ELISA allows for the precise measurement of endogenous levels of natural and/or recombinant CCL21 proteins within the range of 63 pg/mL-4000 pg/mL.