



User's Manual

Mouse CCL6 ELISA Kit



DEIA-XY105



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

Enzyme-linked immunosorbent assay (ELISA) for the determination of mouse CCL6 in cell culture supernatants, serum, or plasma.

Reagents And Materials Provided

1. 96-well plate coated with anti-Mouse CCL6 antibody: 12 x 8 strips
2. Mouse CCL6 standard: 2 x 10 ng
3. Biotin-labeled detection antibody (100x): 130 µL
4. Streptavidin-HRP (100x):130 µL
5. Standard/Sample diluent:30 mL
6. Detection antibody diluent: 12 mL
7. Streptavidin-HRP diluent: 12 mL
8. Wash buffer (20x): 30 mL
9. TMB substrate solution: 12 mL
10. Stop solution: 12 mL
11. Plate adhesive strips: 3 Strips
12. Technical manual: 1 Manual

Storage

2-8°C

Detection Range

15.6 pg/mL-1000 pg/mL

Sensitivity

The minimum detectable dose of mouse CCL6 is typically less than 10 pg/mL.

Specificity

The mouse CCL6 ELISA allows for the detection and quantification of endogenous levels of natural and/or recombinant mouse CCL6 proteins within the range of 15.6 pg/mL-1000 pg/mL.