



User's Manual

Human AIMP1 ELISA Kit

REF

DEIA-XY106



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

Enzyme immunoassay (ELISA) for the determination of Aminoacyl tRNA synthetase complex interacting multifunctional protein 1 (AIMP1) in human serum and cell lysate.

Storage

4°C, please see protocol for component specific storage instructions.

Precision

Intra-Assay Precision: 3.36%-4.5%

Inter-Assay Precision: 12.31%-18.65%

Detection Range

0.625-40 ng/mL

Sensitivity

0.29 ng/mL

Specificity

This kit exhibits no detectable cross-reactivity with GRS, HRS, KRS, NRS, SRS, WRS and YRS.