



User's Manual

Rat Angiotensinogen (Total) ELISA Kit



DEIA-XY109



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

Solid phase sandwich ELISA for the determination of rat Angiotensinogen in serum, EDTA-plasma, urine or cell culture media.

Reagents And Materials Provided

1. Precoated plate: 1 plate
2. Labeled antibody: 1 x 0.4 mL
3. Standard: 2 x 1.0 mL
4. EIA buffer: 1 x 30 mL
5. Solution for labeled antibody: 1 x 12 mL
6. Chromogen: 1 x 15 mL
7. Stop solution: 1 x 12 mL
8. Wash buffer (40x): 1 x 50 mL

Storage

2-8°C

Precision

Intra-Assay Precision: 3.3%-10.0%

Inter-Assay Precision: 8.5%-15.8%

Detection Range

0.08-5 ng/mL

Sensitivity

0.01 ng/mL



Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)



Fax: 1-631-938-8221



Email: info@creative-diagnostics.com