



User's Manual

Human ANGPTL2 ELISA Kit



DEIA-XY110



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

Solid phase sandwich ELISA for the determination of human ANGPTL2 in serum, EDTA-plasma and cell culture supernatant.

Storage

2-8°C

Precision

Intra-Assay Precision: 3.9%-5.9%

Inter-Assay Precision: 6.3%-10.5%

Detection Range

0.05-3.5 ng/mL

Sensitivity

0.01 ng/mL