



User's Manual

Mouse c-MPL/TPOR ELISA Kit



DEIA-XY112



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

Solid phase sandwich ELISA for the determination mouse soluble c-MPL/TPOR in serum and cell culture supernatant.

Reagents And Materials Provided

1. Precoated plate: 1 plate
2. Labeled antibody: 1 x 0.4 mL
3. Standard: 2 x 0.5 mL
4. EIA buffer: 1 x 30 mL
5. Solution for labeled antibody: 1 x 12 mL
6. Chromogen: 1 x 15 mL
7. Stop solution: 1 x 12 mL
8. Wash buffer (40x): 1 x 50 mL

Storage

2-8°C

Precision

Intra-Assay Precision: 4.4%-6.3%

Inter-Assay Precision: 2.9%-6.3%

Detection Range

21.88-1400.0 pg/mL

Sensitivity

3.5 pg/mL



Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)



Fax: 1-631-938-8221



Email: info@creative-diagnostics.com