



## User's Manual

# Coxsackievirus IgM ELISA Kit



DEIA-XY115



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

Enzyme immunoassay (ELISA) for the determination of human IgM antibodies directed against Coxsackie viruses in serum or plasma (EDTA, citrate, heparin).

### Storage

2-8°C

### Precision

Intra-Assay Precision: 2.4%-7.5%

Inter-Assay Precision: 3.0%-12.1%

### Sensitivity

>90%

### Specificity

>90%



Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)



Fax: 1-631-938-8221



Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)