



User's Manual

Legionella pneumophila 1-7 IgG ELISA Kit



DEIA-XY124



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

Enzyme immunoassay (ELISA) for the determination of human IgG antibodies directed against Legionella pneumophila serotypes 1 to 7 in serum or plasma (EDTA, citrate, heparin).

Storage

2-8°C

Precision

Intra-Assay Precision: 4.5%-5.0%

Inter-Assay Precision: 5.4%-7.3%