



User's Manual

Human deamidated gliadin protein IgA/IgG ELISA Kit



DEIA-XY130



96T





This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

Immunometric enzyme immunoassay (ELISA) for the determination of antibodies (IgA and IgG) against deamidated gliadin protein epitopes (DGP) in human serum or plasma.

Storage

2-8°C

Precision

Intra-Assay Precision: 3.2%-5.8%

Inter-Assay Precision: 2.1%-4.9%

Detection Range

6.3-100 U/mL



Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)



Fax: 1-631-938-8221



Email: info@creative-diagnostics.com