



User's Manual

6-Keto-PGF1-alpha ELISA kit



DEIA-XY2035



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

This is an ELISA for the quantitative analysis of 6-k-PGF1a levels in biological fluid.

Storage

Store kit at 4°C at all times.

Specificity

6- keto-Prostaglandin F1a: 100.0%

Prostaglandin F1a: 76.77%

2, 3-dinor-6-keto-Prostaglandin F1a: 64.48%

6- keto-Prostaglandin E1: 49.92%

Prostaglandin F2a: 29.94%

Prostaglandin E2: 11.84%

Prostaglandin D2: 8.36%

11- d-Thromboxane B2: 1.07%

15- k-Prostaglandin F2a: 0.66%

13, 14-d-15-k-prostaglandin F2a: 0.48%

Prostaglandin A2: 0.30%

Prostaglandin A1: 0.25%

Thromboxane B2: 0.15%

Prostaglandin B2: 0.02%

Prostaglandin B1: 0.02%

Leukotriene B4: 0.01%