



## User's Manual

# Cyclic AMP ELISA kit



DEIA-XY2038



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

This is a competitive ELISA for the quantitative analysis of cAMP levels in biological fluid.

### Storage

Store lyophilized conjugate at -20°C or less.

Store all other kit components at 4°C at all times.

### Specificity

Cyclic AMP: 100.0%

Cyclic GMP: 0.07%

Adenosine: <0.01%

Adenine: <0.01%

AMP: <0.01%

ADP: <0.01%

ATP: <0.01%

Guanine: <0.01%

GMP: <0.01%

GTP: <0.01%