



User's Manual

# Rat Prorenin/Renin Total Antigen ELISA kit

REF

DEIA-XY2040



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

This is an ELISA (Enzyme-Linked Immunosorbent Assay) for the quantitative analysis of total rat prorenin/renin in biological fluid.

### Storage

1. Store this kit and its components at 4°C until use.
2. Reconstituted standards and primary may be stored at -70°C for later use. DO NOT freeze/thaw the Standards and Primary Antibody more than once.

### Detection Range

Typical rat prorenin levels range from 0 to 400 ng/mL. The assay range for this ELISA is 0.1 to 100 ng/mL. Samples with prorenin levels above 100 ng/mL should be diluted in Blocking Buffer and retested.