



User's Manual

Rat tPA ELISA Active ELISA kit

REF

DEIA-XY2041



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

This is an ELISA (Enzyme-Linked Immunosorbent Assay) for the quantitative analysis of active tPA levels in biological fluid.

Storage

1. Store this kit and its components at 4°C until use.
2. Do not freeze.

Detection Range

0.05 - 10 ng/mL