



User's Manual

Human Factor XI Total Antigen ELISA kit

REF

DEIA-XY2044



96T



RUO



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

This assay is intended for the quantitative determination of total Factor XI antigen in human plasma.

Storage

1. Store this kit and its components at 4°C until use.
2. The reconstituted Standard and Primary Antibody may be stored at -70°C for later use. DO NOT freeze/thaw the Standard or Primary Antibody more than once.

Detection Range

3.0 - 6.0 µg/mL