



**User's Manual**

# Human Elafin (Trappin-2) ELISA kit

**REF**

**DEIA-XY2150**



**96T**

**RUO**

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

**Creative Diagnostics**

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)**  **Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)**

---

## PRODUCT INFORMATION

### Intended Use

The human elafin ELISA has been developed for the quantitative measurement of natural and recombinant elafin in serum, plasma and cell culture medium. In serum or plasma samples elafin can be measured accurately if serum or plasma samples are diluted at least 10 times.

### Storage

Product should be stored at 4°C. Under recommended storage conditions, product is stable for one year.

### Detection Range

878 - 10,000 pg/mL

### Detection Limit

878 pg/mL

