



## User's Manual

# Human Ficolin-1 ELISA kit



DEIA-XY2155



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

The human Ficolin-1 ELISA kit is to be used for the in vitro quantitative determination of human Ficolin-1 in serum and plasma samples. This kit is intended for laboratory research use only and is not for use in diagnostic or therapeutic procedures. The analysis should be performed by trained laboratory professionals.

### Storage

Product should be stored at 4°C. Under recommended storage conditions, product is stable for one year.

### Detection Range

3.1 - 200 ng/mL

### Detection Limit

3.1 ng/mL