



## User's Manual

# Human Ficolin-3 ELISA kit

REF

DEIA-XY2156



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

The human Ficolin-3 ELISA kit is to be used for the in vitro quantitative determination of human Ficolin-3 in serum, plasma and cell culture supernatant samples. This kit is intended for laboratory research use only and is not for use in diagnostic or therapeutic procedures. The analysis should be performed by trained laboratory professionals.

### Storage

Product should be stored at 4°C. Under recommended storage conditions, product is stable for one year.

### Detection Range

3 - 500 ng/mL

### Detection Limit

3 ng/mL

