



User's Manual

Mouse S100A9/Calprotectin ELISA kit



DEIA-XY2170



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The mouse S100A9 ELISA kit is to be used for the in vitro quantitative determination of mouse S100A9 in serum and plasma samples. The mouse S100A9 ELISA will detect murine S100A9 homodimers as well as higher molecular oligomers, including S100A8/A9 multimers.

Storage

Product should be stored at 4°C. Under recommended storage conditions, product is stable for one year.

Detection Range

0.16 - 10 ng/mL

Detection Limit

0.16 ng/mL