



User's Manual

Human sTREM-1 ELISA kit



DEIA-XY2174



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The human sTREM-1 ELISA kit is to be used for the in vitro quantitative determination of human sTREM-1 in serum and plasma samples.

Storage

Product should be stored at 4°C. Under recommended storage conditions, product is stable for one year.

Precision

Intra Assay Precision CV%: < 20%

Inter Assay Precision CV%: < 25%

Detection Range

31 - 2000 pg/mL

Detection Limit

31 pg/mL