



User's Manual

Human TRACP5a ELISA kit



DEIA-XY2176



96T




This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The human TRACP5a ELISA kit is to be used for the in vitro quantitative determination of human TRACP5a in plasma and serum samples.

Storage

Product should be stored at 4°C. Under recommended storage conditions, product is stable for one year.

Detection Range

78 - 5000 pg/mL

Detection Limit

78 pg/mL