



## User's Manual

# Rat MPO ELISA kit



DEIA-XY2177



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

The rat MPO ELISA has been developed for the quantitative measurement of natural and recombinant rat MPO in cell culture medium, plasma and tissue homogenates. In plasma samples MPO can be measured accurately if plasma samples are diluted at least 4 times.

### Storage

Product should be stored at 4°C. Under recommended storage conditions, product is stable for one year.

### Detection Range

3.9 - 250 ng/mL

### Detection Limit

3.9 ng/mL