



User's Manual

Human ADAM15 ELISA Kit



DEIA-XY2180



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

This Human ADAM15 ELISA Kit contains the necessary components required for the quantitative measurement of natural and recombinant human ADAM15 from serum and plasma in a sandwich ELISA format.

Storage

2-8°C

Precision

Intra-assay Precision: 6-8%; Inter-assay Precision: 8-12%

Detection Range

47-3000 pg/ml

Sensitivity

10 pg/mL

Specificity

Human ADAM15

