



## User's Manual

# Human Total Adiponectin ELISA Kit

REF

DEIA-XY2181



96T

RUO



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

Human Adiponectin (Total) immunoassay is a solid phase ELISA designed to measure total (low, middle, and high molecular weight) human Adiponectin in cell culture supernates, serum, and EDTA plasma.

### Storage

2-8°C

### Precision

Intra-assay Precision: 4-8%; Inter-assay Precision: 8-12%

### Detection Range

62.5-4000 pg/ml

### Sensitivity

31 pg/ml

### Specificity

Human ADIPONECTIN only