



User's Manual

Human Soluble Anterior Gradient 3 ELISA Kit

REF

DEIA-XY2183



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

This Human soluble Anterior Gradient 3 (AGR3)/Breast cancer membrane protein 11 ELISA Kit contains the necessary components required for the quantitative measurement of recombinant and/or natural human AGR3 from cell culture supernates, serum and plasma in a sandwich ELISA format.

Storage

2-8°C

Precision

Intra-assay Precision: 4-6%; Inter-assay Precision: 8-10%

Detection Range

39-2500 pg/mL

Sensitivity

4 pg/mL

Specificity

Human sAGR3