



User's Manual

Human Angiopoietin-4 ELISA Kit

REF

DEIA-XY2184



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

Human Angiopoietin 4 immunoassay is a solid phase ELISA designed to measure human Angiopoietin 4 in cell culture supernates, EDTA plasma and serum.

Storage

2-8°C

Precision

Intra-assay Precision: 4-6%; Inter-assay Precision: 8-12%

Detection Range

125-8000 pg/mL

Sensitivity

30 pg/mL

Specificity

Human Angiopoietin 4