



User's Manual

Human Cadherin-16/KSP-Cadherin ELISA Kit



DEIA-XY2187



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

This Human Cadherin-16 (CDH-16) ELISA Kit contains the necessary components required for the quantitative measurement of natural and recombinant human CDH-16 from serum and plasma in a sandwich ELISA format.

Storage

2-8°C

Precision

Intra-assay Precision: 6-8%; Inter-assay Precision: 8-12%

Detection Range

109-7000 pg/ml

Sensitivity

45 pg/mL

Specificity

Human CDH-16