



User's Manual

Human Chromogranin A (19-131) ELISA Kit



DEIA-XY2188



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

For the quantitative determination of vasostatin-2 concentrations in human serum and EDTA plasma.

Storage

2-8°C

Precision

Intra-assay Precision: 6-8%; Inter-assay Precision: 12-14%

Detection Range

0.64-400 ng/mL

Sensitivity

0.32 ng/mL

Specificity

Human Vasostatin-2