



User's Manual

Human Total Betatrophin ELISA Kit

REF

DEIA-XY2191



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

This Human Betatrophin (Total) ELISA Kit contains the necessary components required for the quantitative measurement of recombinant and/or natural human the C-terminal Betatrophin (100-198) and full length betatrophin from serum and plasma in a sandwich ELISA format.

Storage

2-8°C

Precision

Intra-assay Precision: 4-6%; Inter-assay Precision: 8-10%

Detection Range

781-50000 pg/mL

Sensitivity

120 pg/mL

Specificity

Human the C-terminal Betatrophin (100-198) and Betatrophin (21-198)