



User's Manual

Mouse/Rat/Human CTRP10 ELISA Kit

REF

DEIA-XY2192



96T



RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

Mouse/Rat/Human CTRP10 immunoassay is a 3.5 - 4.5 hour solid phase ELISA designed to measure Mouse/Rat/Human CTRP10 in serum and plasma.

Storage

2-8°C

Precision

Intra-assay Precision: 6-8%; Inter-assay Precision: 8-12%

Detection Range

4.88-625 pg/ml

Sensitivity

1 pg/mL

Specificity

Mouse, Rat, Human CTRP10 Globular form

