



User's Manual

Human CTRP12 ELISA Kit



DEIA-XY2193



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

Human CTRP12/FAM132A immunoassay is a 3.5 - 4.5 hour solid phase ELISA designed to measure human CTRP12 in serum and plasma.

Storage

2-8°C

Precision

Intra-assay Precision: 6-8%; Inter-assay Precision: 8-12%

Detection Range

0.049-200 ng/mL

Sensitivity

10 pg/mL

Specificity

Human CTRP12/FAM132A