



User's Manual

Human CTRP7 ELISA Kit



DEIA-XY2201



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

This Human CTRP7 ELISA Kit contains the necessary components required for the quantitative measurement of recombinant and/or natural human CTRP7 from serum and plasma in a sandwich ELISA format.

Storage

2-8°C

Precision

Intra-assay Precision: 6-8%; Inter-assay Precision: 8-12%

Detection Range

625-40000 pg/mL

Sensitivity

100 pg/mL

Specificity

Human CTRP7

