



User's Manual

Mouse DLL4 ELISA Kit



DEIA-XY2204



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

This Mouse Delta-Like Protein 4 (DLL4) ELISA Kit contains the necessary components required for the quantitative measurement of recombinant and/or natural mouse DLL4 from serum, plasma and cell cultures in a sandwich ELISA format.

Storage

2-8°C

Precision

Intra-assay Precision: 4-6%; Inter-assay Precision: 8-10%

Detection Range

9.3-600 pg/mL

Sensitivity

15 pg/mL

Specificity

Mouse DLL4

