



User's Manual

Human ECM-1 ELISA Kit



DEIA-XY2205



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

This human soluble Extracellular Matrix Protein-1 (ECM-1) ELISA Kit contains the necessary components required for the quantitative measurement of recombinant and/or natural human ECM-1 from serum, plasma and cell cultures in a sandwich ELISA format.

Storage

2-8°C

Precision

Intra-assay Precision: 4-6%; Inter-assay Precision: 8-10%

Detection Range

23.4-1500 pg/mL

Sensitivity

10 pg/mL

Specificity

Human soluble ECM-1