



User's Manual

Human Soluble EphA1 ELISA Kit

REF

DEIA-XY2210



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

This Human soluble Ephrin Type-A Receptor 1 (EphA1) ELISA Kit contains the necessary components required for the quantitative measurement of recombinant and/or natural soluble human EphA1 from cell culture supernates, serum and plasma in a sandwich ELISA format.

Storage

2-8°C

Precision

Intra-assay Precision: 4-6%; Inter-assay Precision: 8-10%

Detection Range

39-2500 pg/mL

Sensitivity

12 pg/mL

Specificity

Human EphA1

