



User's Manual

Mouse Endothelial cell specific molecule-1/Endocan ELISA Kit

REF

DEIA-XY2214



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

For the quantitative determination of mouse esm-1 concentrations in serum and EDTA plasma.

Storage

2-8°C

Precision

Intra-assay Precision: 6-8%; Inter-assay Precision: 10-12%

Detection Range

11.7-750 pg/ml

Sensitivity

2 pg/ml

Specificity

Mouse ESM-1 only