



User's Manual

Human Soluble FAM3B/Pancreatic Derived Factor ELISA Kit

REF

DEIA-XY2215



96T



RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

This Human soluble Pancreatic Derived Factor/FAM3B ELISA Kit contains the necessary components required for the quantitative measurement of recombinant and/or natural human Pancreatic Derived Factor/FAM3B from cell culture supernates, serum and plasma in a sandwich ELISA format.

Storage

2-8°C

Precision

Intra-assay Precision: 4-6%; Inter-assay Precision: 8-10%

Detection Range

6.25-200 ng/mL

Sensitivity

1.5 ng/mL

Specificity

Human FAM3B

