



## User's Manual

# Rat Gas6 ELISA Kit



DEIA-XY2217



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

This Rat Gas6 ELISA Kit contains the necessary components required for the quantitative measurement of recombinant and/or natural rat Gas6 from cell culture supernates, serum and plasma in a sandwich ELISA format.

### Storage

2-8°C

### Precision

Intra-assay Precision: 4-6%; Inter-assay Precision: 8-10%

### Detection Range

39-2500 pg/mL

### Sensitivity

10 pg/mL