



User's Manual

Human Soluble JAM-A ELISA Kit

REF

DEIA-XY2219



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address:** 45-1 Ramsey Road, Shirley, NY 11967, USA

 **Tel:** 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  **Fax:** 1-631-938-8221

 **Email:** info@creative-diagnostics.com  **Web:** www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

For the quantitative determination of human soluble jam-a/cd321 concentrations in serum and plasma.

Storage

2-8°C

Precision

Intra-assay Precision: 4-6%; Inter-assay Precision: 8-12%

Detection Range

7.8-500 pg/ml

Sensitivity

3 pg/mL

Specificity

Human soluble JAM-A