



User's Manual

Human Kallikrein 13 ELISA Kit



DEIA-XY2220



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

This Human Kallikrein 13/KLK13 ELISA Kit contains the necessary components required for the quantitative measurement of recombinant and/or natural human kallikrein 13 from serum and in a sandwich ELISA format.

Storage

2-8°C

Precision

Intra-assay Precision: 4-6%; Inter-assay Precision: 8-10%

Detection Range

39-2500 pg/mL

Sensitivity

10 pg/mL

