



User's Manual

# Mouse/Rat Lipocalin-13 ELISA Kit



DEIA-XY2223



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

This Mouse/Rat Lipocalin 13 ELISA Kit contains the necessary components required for the quantitative measurement of recombinant and/or mouse lipocalin 13, serum and plasma in a sandwich ELISA format.

### Storage

2-8°C

### Precision

Intra-assay Precision: 4-6%; Inter-assay Precision: 8-10%

### Detection Range

78-5000 pg/mL

### Sensitivity

30 pg/mL

### Specificity

Mouse, Rat

